

The World's 1st A.I.-Equipped • Portable Artificial Nose (A.I. Nose)

~ Portable Odor Detector/Identifier] “Aino-Pro” Development Press Release~

※1 (Current Research as of January 15, 2018)

The Shizuoka Cancer Center Research Institute (Shizuoka Prefecture Sunto-gun Nagaizumi-cho, President: Ken Yamaguchi and Deputy Director: Masatoshi Kusuhara) and Tokai-Denshi Inc. (Headquarters: Fuji City, representative: Kazunari Sugimoto) jointly developed an A.I.-based Portable Odor detector/Identifier which analyzes malodorous components emitted from cancer patients.

## 1. Background

### Issues in the Medical Field

During the advanced stages of cancer, especially in patients with uterine and breast cancer, occasionally there is a strong odor or “sick smell”. Patients often can become annoyed by the odor, suffer pain, and often times their medical care life is extremely difficult.

If there is a device that can analyze the disease odor elements easily with high accuracy, it can help doctors select the optimum deodorizing method (the method of treating the affected part and the type of deodorant). This will allow for specialists to choose which treatments and how best to apply them to the care of each patient.



Photo 1: Shizuoka Cancer Center Research Institute



Photo 2: Aino-Pro Exterior

Therefore, this device is made possible by the cooperation of the Pharma Valley Center (Director: Katsunori Ueda) at the Ministry of Economy, Trade and Industry for the FY 2013 and 2014 medical device development project (comprehensive special zoning promotion commission fee) "for odor measurement using odor sensor", The Shizuoka Cancer Center Research Institute (deputy director Masatoshi Kusuhara), and Tokai-Denshi Inc.

Up until now, odor measuring devices and the like have been on the market. However, since there was no device capable of identifying the kind of odor like human olfaction, until now patients have had to rely on the smell of doctors and nurses. At this time, the Portable Odor Detector/Identifier developed by Tokai-Denshi Inc. has made it possible to identify dozens of odors. The device consists of artificial intelligence and a portable simple gas chromatograph, and the A.I. inside the device learns the odor component, thereby making it possible to distinguish odors.

A.I. inside the device learns the odor component, thereby making it possible to distinguish odors.

## 2. Summary

People can easily distinguish odors drifting through the air. This is largely due in part to the work of about 400 olfactory cells. Odor is absorbed by olfactory cells and is replaced with electric signals and transmitted to the brain. This time, the disease odor measuring device began with a survey of gas sensors that can distinguish odors and VOC gases (volatile organic compounds). However, the biggest challenge was it is difficult to obtain the sensitivity of human olfactory cells and the ability to discriminate odors. Therefore, we developed a highly sensitive and portable gas analyzer (simple gas chromatograph) capable of separating and identifying nearly 100 kinds of VOC gases, and solved the problem by installing artificial intelligence.

Through using artificial intelligence, we discover individual features and weak changes. This makes it possible to identify odors by learning the characteristics of odors in each. As a sensor that detects odors, a high sensitivity sensor that can detect 5 ppb to 200 ppm using a non-concentrated gas was used. This was using a PID sensor (photoionization detector) and a semiconductor gas sensor. In addition, by using a 30 m long gas analysis tube (capillary column), the odor gas was separated and the analytical performance was improved. The odor is analyzed by searching for features using artificial intelligence detection signals obtained from the surrounding air.




Photo 3 : Aino-Pro System Configuration

The device starts with “learning” the sick odor as well as mimicking the human olfactory sense in order to perform accurate odor determination. We acknowledge that odor detection can be used not only in the medical field but also in various other situations. As such, we will develop devices specifically tailored to categorical odor identification rather than universal odor identification. By studying by category, since the kinds of gases that specify odors are almost the same, it is expected that the discrimination ability (correct answer rate) will be further improved. In other words, odor classification category detection can be customized to the device such as offensive odor from a livestock farm, a poultry farm, a pig farm, a food processing factory, air pollution of chemical plant, scent development such as food sample, halitosis measurement, drug or doping inspection etc.

Looking forward to the future, we are looking to prepare a free monitoring machine and recruit collaborative research laboratories or companies while we continue to develop high sensitivity judgment and super sensitivity sensors exceeding human olfactory ability. By 2020 we plan to complete the next

generation artificially intelligent electronic nose for commercialization in the medical, drug, pollution, food, and aroma industry.

This equipment will be exhibited February 2<sup>nd</sup>, 2018 at "Mt. Fuji Industry-Academia-Government-Finance Collaboration Forum 2018" at Plaza Verde, 1-1-4 Otemachi, Numazu City.





Photo 4 : Aino-Pro Internal photo during an experiment

Photo 5 : (from left) M.Kusuhara,M.D., Ph.D. and Developers

### 3. Product Specs

#### Aino-Pro

Gas Sensor	PID (Photoionization detector) ,Semiconductor gas sensor Sensitivity 5 ppb~200 ppm (Isobutylene)
Collecting Pipe Measurement	Can measure approximately up to 200 times concentration measurable up to 25 ppt (Isobutylene)
Carrier Gas	Indoor air
Sample Gas Measurement	Available for Sample bag, Surrounding gas, Thermal desorption tube
Measurement Sequence	Up to 19 different measurement sequences can be set (Column Flow Rate, Makeup Flow Rate, Column Temperature Program, Thermal Desorption Tube Temperature Program)
Backflush Function	Built-in 8-way valve utilizing a backflush method to avoid column contamination
A.I.	It can be set by category of odor (number of input variable nodes, number of hidden layer nodes, learning rate, activation function, parameter update, overfitting control)
Judgement Output	A user can make this system learn odor type through 50 ways or 6 stages of odor intensity.
Easy Operation	It has 20 kinds of measurement sequences and can be selected on the PC screen or measured with the dial type measurement menu..
PC-less Measurement	In the case of PC-less measurement, after downloading from a USB, PC analysis is possible.
Time Table Measurement	With the time table function, you can program measurements of up to 9 types of

	daily, hourly, constant time interval etc.
Consumable Goods※1	VOC filter, capillary column, various heaters, and pumps
Size and Weight	210 w × 450 d × 320 h, Weight Approximately 10 kg
Low Power Consumption	100 V 350 W or less (Standby time at 20 W, Operating time at 300 W, Maximum 350 W)

#### 4. System Configuration

- **Aino-Pro** ( Portable Odor Detector/Identifier)
- A.I. analysis software standard
- Windows PC for data collection analysis (Windows 10)...User arrangement

#### 5. Practical Use

- Odor analysis
- Environmental analysis
- Halitosis analysis

#### 6. Release Date Complete order production system

Orders will begin on March 1, 2018, shipping will begin from July (scheduled)

Selling Price 3.9 million yen (planned) standard package will include the **Aino-Pro** main body and A.I. analysis software

A separate maintenance contract is required.

(It depends on usage form, existence of version up contract etc.)

(※ 1. Consumables , fees, and costs will be charged separately)

---

★★For inquiries or requests for more information★★

Tokai Electronics Co., Ltd. Medical System Equipment Division 2-34-13 Akebonocho, Tachikawa-shi,

Tokyo Olympic Building No. 3 Room 203 TEL:042-526-0905/FAX:042-526-0906

E-mail: [info@tokai-denshi.co.jp](mailto:info@tokai-denshi.co.jp) URL : [www.tokai-denshi.co.jp](http://www.tokai-denshi.co.jp)

---